



RETRACTOR MAKER

PRODUCT BROCHURE OF

Harvey Jackson Self Retaining Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Harvey Jackson Self Retaining Retractor

SKU RM-HJSRR-1

\$229 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
Harvey Jackson Self Retaining Retractor	28 cm overall length, 45 × 32 mm jaws, 3×3 sharp prongs — 1	\$229
Protective Sleeve / Box	Individual transit packaging — 1	\$229

PRODUCT DESCRIPTION

The Harvey Jackson Self Retaining Retractor is a heavy-duty surgical retractor engineered specifically to provide powerful, hands-free soft-tissue and muscular retraction during deep spinal, orthopedic, and pelvic operations. Designed with a robust mechanical frame, this instrument maintains broad, unyielding exposure across thick muscular planes, eliminating manual assistant fatigue during extensive bone decompressions.

Built to provide optimal surface area exposure in deep cavities, this Harvey Jackson Self Retaining Retractor configuration features an overall tool length of 28 cm (\$229). Its working jaws measure 45 mm in depth by 32 mm in width, incorporating a 3×3 outward-curving sharp prong layout that anchors securely into dense fibrous fascia and muscle. The integrated multi-position mechanical ratchet lock holds targeted tissue counter-tension automatically, preventing accidental slippage throughout lengthy surgical procedures.

Forged from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with a reflection-free satin matte finish, Retractor Maker manufactures the Harvey Jackson Self Retaining Retractor to withstand heavy mechanical loads and repeated hospital autoclave cycles, providing spinal, orthopedic, and neurosurgical operating theaters worldwide with exceptional durability and clinical stability.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)