



RETRACTOR MAKER

PRODUCT BROCHURE OF

Holzhimer Self-Retaining Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Holzimer Self-Retaining Retractor

SKU RM-HSRR-4IN-1

\$49 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
Holzimer Self-Retaining Retractor — 4" (10 cm) length, 4×4 outward prongs (Sharp or Blunt)	RM-HSRR-4IN-1-1	\$49

PRODUCT DESCRIPTION

The Holzimer Self-Retaining Retractor is an indispensable minor surgery retractor engineered specifically to provide dependable, hands-free wound edge separation and stable tissue exposure during superficial surgical procedures. Designed with an ultra-compact chassis, this retractor allows surgeons to anchor wound margins securely, freeing assistant hands while preserving an unobstructed line of sight in tight operative spaces.

Built in a unified space-saving profile to fit confined operating arenas, this Holzimer Self-Retaining Retractor model features an overall instrument length of 4 inches (10 cm, \$49). Its working ends incorporate a 4×4 outward-curving prong layout designed to hook firmly behind cutaneous borders and superficial fascial layers. Surgeons can customize the instrument with selectable Sharp tips for firm fibrous tissue purchase or Blunt tip variants to safeguard delicate neurovascular structures. The integrated precision thumb-screw drive ensures smooth, micro-adjustable arm spreading that stays locked under counter-traction without accidental slipping.

Forged from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with a non-reflective satin matte finish, Retractor Maker manufactures the Holzimer Self-Retaining Retractor to withstand repeated hospital autoclaving cycles, providing plastic, hand, and general surgery clinics with unyielding structural stability, corrosion resistance, and clinical longevity.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)