



RETRACTOR MAKER

PRODUCT BROCHURE OF

Breisky Vaginal Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.

Donald Cameron



Breisky Vaginal Retractor

SKU RM-BVR-1

\$89 USD

COMPONENT LIST

QUANTITY	COMPONENT	DIMENSIONS / DESCRIPTION
1	100 mm × 25 mm — 28 cm to 36 cm	Superficial vaginal and introitus exposure — \$89.00
1	100 mm × 35 mm / 40 mm — 28 cm to 36 cm	Broad anterior/posterior vaginal wall holding — \$89.00
1	130 mm × 20 mm / 30 mm — 28 cm to 36 cm	Intermediate vaginal vault exposure — \$94.00
1	130 mm × 40 mm — 28 cm to 36 cm	Deep pelvic floor reconstruction — \$94.00
1	140 mm × 35 mm — 28 cm to 36 cm	Vaginal hysterectomy cervical access — \$99.00
1	160 mm × 20 mm / 35 mm — 28 cm to 36 cm	Deep vault access in parous anatomy — \$99.00
1	180 mm × 25 mm — 28 cm to 36 cm	Maximum-depth pelvic and apical access — \$104.00
1	Protective Sleeve / Box	Standard packaging — Included

PRODUCT DESCRIPTION

The Breisky Vaginal Retractor is a premier gynecological and pelvic reconstructive surgical retractor engineered specifically to provide deep, unobstructed visualization of the vaginal vault and cervix. Designed with an elongated, smooth spatula blade that curves gracefully at its distal tip, this instrument depresses the vaginal wall and displaces surrounding soft tissues without pinching or lacerating mucosa.

Built to accommodate variable patient anatomy and deep pelvic procedures, the Breisky Vaginal Retractor is offered in selectable overall lengths of 28 cm (\$89), 31 cm (\$94), 33.5 cm (\$99), and 36 cm (\$104), with extended production runs accommodating profiles up to 15.25 inches (38.73 cm). Surgeons can choose from 10 distinct blade dimensions spanning 100 × 25 mm to 180 × 25 mm. Its ergonomically angled, flat grooved handle provides surgical assistants with steady counter-traction and non-slip control during long pelvic operations.

Forged from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with a reflection-free satin matte finish, Retractor Maker manufactures the Breisky Vaginal Retractor to withstand repetitive hospital sterilization cycles, providing obstetrics, gynecology, and urogynecology teams with high corrosion resistance and dependable durability.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)