



RETRACTOR MAKER

PRODUCT BROCHURE OF

Tessier Malleable Ribbon Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Tessier Malleable Ribbon Retractor

SKU RM-TMRR-22-1

\$44 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
Narrow Blade — 22 cm length, narrow malleable spatula — Orbital floor exploration, medial canthal access, pediatric craniofacial	RM-TMRR-22-1-1	\$44
Medium Blade — 22 cm length, medium malleable spatula — Standard orbital rim exposure, zygomatic arch and maxillary osteotomies	RM-TMRR-22-1-2	\$44
Wide Blade — 22 cm length, wide malleable spatula — Broad forehead flap elevation, coronal approach, and globe shielding	RM-TMRR-22-1-3	\$44
Set of Three Package — Combined Narrow, Medium & Wide — Complete craniofacial malleable tray set for all facial bone corridors	RM-TMRR-22-1-4	\$94

PRODUCT DESCRIPTION

The Tessier Malleable Ribbon Retractor is a specialized craniofacial, plastic, and orbital surgical retractor engineered specifically to provide safe, non-traumatic soft-tissue displacement and globe protection during complex facial reconstructions. Designed with a ductile, malleable ribbon chassis, this instrument allows surgeons to contour the blade intraoperatively around the orbital rim, zygoma, and maxilla, providing stable tissue retraction without causing focal pressure necrosis on delicate periocular structures.

Built to withstand frequent autoclave sterilization runs while remaining highly resistant to corrosion, this Tessier Malleable Ribbon Retractor configuration features an overall tool length of 22 cm (\$44). Handcrafted from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic), the instrument incorporates a base blade diameter/thickness of 9 mm with a physical instrument weight of 0.100 kg that provides optimal tactile sensitivity. The line is highly customizable across three distinct width setups—Narrow Blade layout, Medium Blade layout, and Wide Blade layout—or selectable as a comprehensive Set of Three package upgrade combining all three configurations for \$94 (+\$50.00 option upgrade).

Forged with an anti-glare satin matte finish, Retractor Maker manufactures the Tessier Malleable Ribbon Retractor to withstand repeated intraoperative bending and hospital reprocessing cycles, providing maxillofacial suites, plastic surgery clinics, and craniofacial centers worldwide with exceptional ductility, corrosion resistance, and complete autoclave durability.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)