



RETRACTOR MAKER

PRODUCT BROCHURE OF

Femoral Calcar Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Femoral Calcar Retractor

SKU RM-FCR-1

\$44 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
Femoral Calcar Retractor	25 mm tip width, 5.08 cm (2") depth, ring handle	\$44
Protective Sleeve / Box	Individual transit packaging	\$44

PRODUCT DESCRIPTION

The Femoral Calcar Retractor is an indispensable orthopedic surgical instrument engineered specifically to provide clear, stable exposure of the medial femoral neck and calcar during total hip arthroplasty and femoral head resection. Designed to hook over the calcar femorale, this instrument allows surgical assistants to retract the medial capsular flap and surrounding adductor musculature safely away from femoral broaches.

Built for optimal joint space exposure, this Femoral Calcar Retractor configuration (\$44) features a standardized 25 mm blade tip width paired with a 5.08 cm (2-inch) blade depth. The curved distal tip is contoured to wrap smoothly around the calcar bone cortex, providing exceptional counter-traction stability without slipping into the obturator vascular zone. Its handle incorporates an ergonomic smooth ring end that provides a secure manual grip and reduces hand fatigue during continuous femoral preparation.

Forged from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with a reflection-free satin matte finish, Retractor Maker manufactures the Femoral Calcar Retractor to withstand heavy mechanical torque and repeated hospital autoclave cycles, providing orthopedic joint replacement surgeons worldwide with unyielding clinical durability and precision.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractor maker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractor maker.com
retractor manufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)