



RETRACTOR MAKER

PRODUCT BROCHURE OF

Bailey Rib Contractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Bailey Rib Contractor

SKU RM-BRC-1

\$159 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
Bailey Rib Contractor	20 cm overall length, claw-tip mechanical compressor	\$159
Protective Sleeve / Box	Individual transit packaging	\$159

PRODUCT DESCRIPTION

The Bailey Rib Contractor (also known as the Bailey-Gibbon rib approximator) is a specialized thoracic surgical instrument engineered to approximate adjacent ribs following a thoracotomy. By bringing spread ribs together under controlled mechanical compression, this instrument allows surgical teams to place and tighten pericostal sutures without tension or bone fracture.

Designed for robust anatomical engagement, this Bailey Rib Contractor model features an overall length of 20 cm, providing adequate reach across thoracotomy closure sites. The working ends incorporate strong, curved claw-like tips specifically contoured to seat securely over the superior and inferior borders of adjacent rib segments without slipping. Its mechanical thumbscrew mechanism allows the surgeon to finely tune compressive force, holding the chest wall in anatomical reduction during closure.

Forged from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with a reflection-free satin matte finish, Retractor Maker manufactures the Bailey Rib Contractor to withstand high compressive torque and repeated steam sterilization cycles, providing cardiothoracic and trauma teams with long-lasting clinical dependability.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)