



RETRACTOR MAKER

PRODUCT BROCHURE OF

Fiber optic Deaver Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Fiber optic Deaver Retractor

SKU RM-FODR-1

\$139 USD

COMPONENT LIST

QUANTITY	COMPONENT	DIMENSIONS / DESCRIPTION
1	Fiber optic Deaver Retractor	8" (20.3 cm) handle, 1" (25 mm) curved blade, fiber-optic bundle
1	Universal Light Cable Adapter	Connects to standard OR fiber-optic cables
1	Protective Sleeve / Box	Standard transit packaging

PRODUCT DESCRIPTION

The Fiber optic Deaver Retractor is a premier illuminated surgical retractor engineered specifically to combine broad, atraumatic visceral retraction with direct, shadowless illumination within deep anatomical cavities. By directing high-intensity light along the inner curve of its blade, this instrument eliminates shadowing from overhead OR lights, allowing surgeons to explore deep pelvic and retroperitoneal spaces with superior optical clarity.

Designed for deep surgical site illumination and cavity clearance, this Fiber optic Deaver Retractor configuration features an 8-inch (20.3 cm) handle (\$139) paired with a 1-inch (25 mm) wide curved blade. The deeply curved blade architecture wraps safely around viscera and muscular walls, providing expansive counter-traction without causing focal organ contusions. Its integrated fiber-optic bundle includes a universal fitting compatible with standard hospital light source cables (Storz, Olympus, Wolf, ACMI).

Forged from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) and high-transmission optical glass fibers with a reflection-free satin matte finish, Retractor Maker manufactures the Fiber optic Deaver Retractor to withstand repetitive hospital autoclave cycles, delivering exceptional durability, optical transmission, and clinical reliability for general, gynecological, and urological surgical suites worldwide.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)