



# RETRACTOR MAKER

## PRODUCT BROCHURE OF

Hohmann Retractor Bent Extra Long Handle

### CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



# Hohmann Retractor Bent Extra Long Handle

**SKU** RM-HRBELH-1

**\$69 USD**

## AVAILABLE SIZES & PRICING

| SIZE / OPTION                             | SKU / DETAILS | PRICE (USD) |
|-------------------------------------------|---------------|-------------|
| 9 ¾" (24.76 cm) — Extra-long bent handle  | RM-HRBELH-1-1 | \$69        |
| 10.0" (25.40 cm) — Extra-long bent handle | RM-HRBELH-1-2 | \$79        |

## PRODUCT DESCRIPTION

The Hohmann Retractor Bent Extra Long Handle is an advanced orthopedic bone leverage instrument engineered specifically to provide superior mechanical advantage and clear visualization during total hip and knee arthroplasty. Designed with an elongated, downward-bent handle stem, this retractor drops the assistant's hands and arms well below the operative plane, ensuring an unobstructed visual line of sight for reaming and component seating.

Built for heavy joint and long bone exposure, the Hohmann Retractor Bent Extra Long Handle features an extra-long handle profile across all versions, allowing assistants to maintain sustained counter-traction with minimal muscle fatigue. Surgeons can select between two primary overall tool lengths: 9 ¾ inches (24.76 cm, \$69) and 10 inches (25.40 cm, \$79). The contoured beak tip hooks securely over cortical bone margins, preventing accidental displacement during high-torque joint manipulation.

Forged from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with a reflection-free satin matte finish, Retractor Maker manufactures the Hohmann Retractor Bent Extra Long Handle to withstand intense mechanical loads and repeated hospital autoclave cycles, delivering exceptional durability, corrosion resistance, and clinical reliability for orthopedic trauma and joint replacement centers worldwide.



# THANK YOU

We look forward to discussing your surgical retractor requirements.

## CONTACT RETRACTOR MAKER

EMAIL

**[info@retractormaker.com](mailto:info@retractormaker.com)**

WHATSAPP

**+92 327 8345101**

WEBSITES

**[retractormaker.com](http://retractormaker.com)**  
**[retractormanufacturer.com](http://retractormanufacturer.com)**

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

## FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)