



RETRACTOR MAKER

PRODUCT BROCHURE OF

Aufricht Nasal Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Aufricht Nasal Retractor

SKU RM-ANR-6IN-1

\$59 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
Without suction (standard)	RM-ANR-6IN-1-1	\$59
With suction tube	RM-ANR-6IN-1-2	\$199

PRODUCT DESCRIPTION

The Aufricht Nasal Retractor is an industry-standard surgical instrument engineered to provide steady elevation of the nasal skin envelope and cartilage during rhinoplasty and reconstructive nasal procedures. Designed with an 8 mm wide flat, solid blade terminating in a curved tip profile, this instrument safely engages deep intranasal tissues while distributing retraction forces evenly across the nasal dorsum.

Engineered for procedural flexibility, this 6-inch (15.2 cm) Aufricht Nasal Retractor is manufactured in two functional configurations: a standard solid blade model (\$59) and an integrated suction tube configuration (\$199). The suction model incorporates a dedicated evacuation channel along the shaft, allowing surgeons to aspirate blood and irrigation fluids from the operative field simultaneously without requiring a separate suction device.

Machined from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with an anti-reflective satin finish, Retractor Maker supplies the Aufricht Nasal Retractor non-sterile and fully autoclavable inside protective packaging, providing plastic, aesthetic, and ENT surgical teams with long-lasting durability and exceptional clinical utility.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractor maker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractor maker.com
retractor manufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)