



# RETRACTOR MAKER

## PRODUCT BROCHURE OF

Image-Line ACF Retractor With Distraction Quick Release  
Set

### CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



# Image-Line ACF Retractor With Distraction Quick Release Set

**SKU** RM-ILACF-QR-1

**\$4,999 USD**

## COMPONENT LIST

| QUANTITY | COMPONENT                           | DIMENSIONS / DESCRIPTION       |
|----------|-------------------------------------|--------------------------------|
| 1        | Depth Gauge                         | Color-coded                    |
| 2        | Spring-loaded Handle                | Quick-Release Mechanism        |
| 1        | Longitudinal Frame                  | Quick-Release Standard         |
| 1        | Longitudinal Frame (Panther Coated) | Glare-Absorbing Surface        |
| 1        | Transverse Frame                    | Quick-Release Standard         |
| 1        | Transverse Frame (Panther Ceramic)  | Ceramic Non-Reflective Coating |
| 2        | Blunt Blade                         | 30 mm, silver                  |
| 2        | Blunt Blade                         | 35 mm, blue                    |
| 2        | Blunt Blade                         | 40 mm, red                     |
| 2        | Blunt Blade                         | 45 mm, green                   |
| 2        | Blunt Blade                         | 50 mm, violet                  |
| 2        | Blunt Blade                         | 55 mm, black                   |
| 2        | Blunt Blade                         | 60 mm, gold                    |
| 2        | Blunt Blade                         | 65 mm, navy                    |

| QUANTITY | COMPONENT                | DIMENSIONS / DESCRIPTION |
|----------|--------------------------|--------------------------|
| 2        | Blunt Blade              | 70 mm, light green       |
| 2        | Blunt Blade              | 75 mm, orange            |
| 1        | Blade, Short Teeth       | 30 mm, silver            |
| 1        | Blade, Short Teeth       | 35 mm, blue              |
| 1        | Blade, Short Teeth       | 40 mm, red               |
| 1        | Blade, Short Teeth       | 45 mm, green             |
| 1        | Blade, Short Teeth       | 50 mm, violet            |
| 1        | Blade, Short Teeth       | 55 mm, black             |
| 1        | Blade, Short Teeth       | 60 mm, gold              |
| 1        | Blade, Short Teeth       | 65 mm, navy              |
| 1        | Blade, Long Teeth        | 30 mm, silver            |
| 1        | Blade, Long Teeth        | 35 mm, blue              |
| 1        | Blade, Long Teeth        | 40 mm, red               |
| 1        | Blade, Long Teeth        | 45 mm, green             |
| 1        | Blade, Long Teeth        | 50 mm, violet            |
| 1        | Blade, Long Teeth        | 55 mm, black             |
| 1        | Blade, Long Teeth        | 60 mm, gold              |
| 1        | Blade, Long Teeth        | 65 mm, navy              |
| 2        | Micro Blade, Blunt       | 14 × 35 mm, blue         |
| 2        | Micro Blade, Blunt       | 14 × 40 mm, red          |
| 2        | Micro Blade, Blunt       | 14 × 45 mm, green        |
| 2        | Micro Blade, Blunt       | 14 × 50 mm, violet       |
| 2        | Micro Blade, Blunt       | 14 × 55 mm, black        |
| 2        | Micro Blade, Blunt       | 14 × 60 mm, gold         |
| 2        | Micro Blade, Short Teeth | 14 × 35 mm               |
| 2        | Micro Blade, Short Teeth | 14 × 40 mm               |

| QUANTITY | COMPONENT                         | DIMENSIONS / DESCRIPTION           |
|----------|-----------------------------------|------------------------------------|
| 2        | Micro Blade, Short Teeth          | 14 × 45 mm                         |
| 2        | Micro Blade, Short Teeth          | 14 × 50 mm                         |
| 2        | Micro Blade, Short Teeth          | 14 × 55 mm                         |
| 2        | Micro Blade, Short Teeth          | 14 × 60 mm                         |
| 1        | CVB Twist Drill                   | 8 mm depth, 5¾" length             |
| 1        | CVB Distractor Drill Guide, Left  | Standard                           |
| 1        | CVB Distractor Drill Guide, Right | Standard                           |
| 1        | Vertebral Body Distractor, Left   | 100 mm rack                        |
| 1        | Vertebral Body Distractor, Right  | 100 mm rack                        |
| 2        | CVB Screw Driver, Silicone Handle | Standard                           |
| 1        | Multilevel Distractor, Left       | 120 mm rack, 90 mm extension       |
| 1        | Multilevel Distractor, Right      | 120 mm rack, 90 mm extension       |
| 1        | Bone Graft Holder                 | 8 × 15 mm tip, 210 mm length       |
| 1        | Bone Graft Impactor               | 3 mm                               |
| 1        | Bone Graft Impactor               | 5 mm                               |
| 1        | Bone Graft Impactor               | 8 mm                               |
| 1        | Graft Cutter                      | Adjustable depth 8.75", 7 mm width |
| 1        | ACF Caliper                       | 9"                                 |
| 1        | Sterilization Box                 | Heavy-Duty Organizer Box           |

## PRODUCT DESCRIPTION

The Image-Line ACF Retractor With Distraction Quick Release Set is an advanced spinal instrument platform engineered to streamline cervical exposure, interbody distraction, and bone graft insertion during anterior cervical fusion (ACDF) procedures. By integrating specialized quick-release mechanical frame locks and spring-loaded handles, this system allows scrub nurses and spine surgeons to swap blades and adjust tension instantly without disrupting the sterile corridor.

Engineered for demanding cervical spine reconstructions, the Image-Line ACF Retractor With Distraction Quick Release Set features four separate frames—including standard and glare-absorbing Panther Coated/Ceramic longitudinal and transverse configurations. The set includes a color-coded depth gauge, 56 retractor blades (blunt tissue, short teeth, long teeth, and micro profiles), left and right 100 mm and 120 mm rack distractors, a

CVB drill/driver system, a bone graft holder, three impactors (3 mm, 5 mm, 8 mm), an adjustable graft cutter, and a 9-inch caliper.

Manufactured from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic), high-durability Panther ceramic coatings, and autoclavable silicone, Retractor Maker supplies the Image-Line ACF Retractor With Distraction Quick Release Set non-sterile and fully autoclavable inside a dedicated heavy-duty sterilization box, providing premium operating room longevity and surgical precision.



# THANK YOU

We look forward to discussing your surgical retractor requirements.

## CONTACT RETRACTOR MAKER

EMAIL

**[info@retractormaker.com](mailto:info@retractormaker.com)**

WHATSAPP

**+92 327 8345101**

WEBSITES

**[retractormaker.com](http://retractormaker.com)**  
**[retractormanufacturer.com](http://retractormanufacturer.com)**

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

## FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)