



RETRACTOR MAKER

PRODUCT BROCHURE OF

Fisch Retractor with Joint

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.

Donald Cameron



Fisch Retractor with Joint

SKU RM-FRWJ-1

\$109 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
2 teeth right / serrated plate left — 17 cm	RM-FRWJ-1-1	\$109
2 teeth left / serrated plate right — 17 cm	RM-FRWJ-1-2	\$109

PRODUCT DESCRIPTION

The Fisch Retractor with Joint is an indispensable otologic and skull base surgical instrument engineered specifically to provide stable, hands-free exposure of the temporal bone, mastoid, and infratemporal fossa. Designed with an articulable mechanical joint, this retractor conforms dynamically to irregular skull contours, eliminating visual obstruction and maintaining a rigid microsurgical corridor for high-speed drilling.

Built for precision in complex anatomical approaches, this Fisch Retractor with Joint configuration features a standard overall tool length of 17 cm (\$109). The working ends incorporate a 3:3 semi-sharp prong layout and are selectable in asymmetrical component configurations: 2 teeth right combined with a serrated plate left, or 2 teeth left combined with a serrated plate right. This asymmetrical arrangement allows the toothed blade to anchor firmly into the retroauricular periosteum while the opposing serrated plate gently retracts the delicate meatal flap without tearing.

Machined from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with a reflection-free satin matte finish, Retractor Maker produces the Fisch Retractor with Joint to withstand repetitive hospital sterilization cycles, providing ENT and lateral skull base surgical suites worldwide with superior mechanical stability, corrosion resistance, and complete autoclave durability.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)