



RETRACTOR MAKER

PRODUCT BROCHURE OF

B-Z Weitlaner Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



B-Z Weitlaner Retractor

SKU RM-BZWR-1

\$339 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
B-Z Weitlaner Retractor	28 cm overall length, 83 mm deep 4×4 blunt prongs, 90° angled	\$339
Protective Sleeve / Box	Individual transit packaging	\$339

PRODUCT DESCRIPTION

The B-Z Weitlaner Retractor is a specialized, deep-access self-retaining surgical retractor engineered to provide powerful, hands-free exposure in demanding orthopedic, spinal, and deep general surgical procedures. Designed with an extended 28 cm overall length, this instrument allows surgical teams to maintain a broad operative corridor in deep cavities without assistant fatigue.

Built for heavy muscular retraction, this B-Z Weitlaner Retractor configuration features a 90-degree angled chassis that keeps the ring handles positioned completely flat against patient contours, freeing the surgical corridor for power tools and microscopes. Its working end is equipped with 83 mm deep blades featuring a 4×4 blunt prong configuration, providing exceptional tissue grip while minimizing trauma to adjacent neural and vascular structures. The integrated mechanical ratchet locks securely at variable widths.

Manufactured from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with a reflection-free satin matte finish, Retractor Maker produces the B-Z Weitlaner Retractor to strict clinical tolerances, ensuring smooth ratchet engagement, high corrosion resistance, and complete autoclave durability.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)