



RETRACTOR MAKER

PRODUCT BROCHURE OF

Pediatric Deaver Retractors

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved




SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.

Donald Cameron



Pediatric Deaver Retractors

SKU RM-PDR-1

\$69 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
17.5 cm (6.9")	RM-PDR-1-1	\$69
20.0 cm (7.9")	RM-PDR-1-2	\$74
22.0 cm (8.7")	RM-PDR-1-3	\$79

PRODUCT DESCRIPTION

Pediatric Deaver Retractors are specialized pediatric surgical instruments engineered specifically to provide safe, atraumatic visceral and abdominal wall displacement during pediatric general, urological, and thoracic procedures. Scaled down to accommodate smaller anatomical cavities, these retractors allow pediatric surgeons to hold back fragile infant liver lobes, spleen, intestines, and pelvic structures without causing focal organ contusions or serosal tears.

Built with a highly adaptive selection matrix, this Pediatric Deaver Retractors series is offered across three overall lengths: 17.5 cm (\$69) for neonates and infants, 20.0 cm (\$74) for pediatric patients, and 22.0 cm (\$79) for adolescent surgical corridors. Handcrafted from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic), surgeons can select from three working widths: 8 mm, 10 mm, and 12 mm. Its curved hook blade profile incorporates smooth, rounded atraumatic borders, while the ergonomic flat handle provides operating assistants with steady manual leverage.

Forged with an anti-glare satin matte finish, Retractor Maker produces Pediatric Deaver Retractors as non-sterile, re-sterilizable operating room devices built to withstand frequent sterilization runs and repetitive autoclave cycles, providing pediatric surgery centers and children's hospitals worldwide with exceptional durability, corrosion resistance, and complete clinical dependability.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractor maker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractor maker.com
retractor manufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)