



RETRACTOR MAKER

PRODUCT BROCHURE OF

Kocher Langenbeck Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Kocher Langenbeck Retractor

SKU RM-KLR-GEN-1

\$44 USD

COMPONENT LIST

QUANTITY	COMPONENT	DIMENSIONS / DESCRIPTION
1	Narrow Blades (e.g., 25×6, 35×8 mm)	Available Lengths & Handles: 20.3 or 22 cm (Ring / Hollow) · Primary Clinical Role: Focused thyroid, tendon, and minor surgical retraction · Price (USD): \$44.00 / \$49.00
1	Medium Blades (e.g., 28×15, 41×14 mm)	Available Lengths & Handles: 20.3 or 22 cm (Ring / Hollow) · Primary Clinical Role: Standard appendectomy, hernia, and orthopedic joint access · Price (USD): \$44.00 / \$49.00
1	Deep Blades (e.g., 70×14, 80×16 mm)	Available Lengths & Handles: 20.3 or 22 cm (Ring / Hollow) · Primary Clinical Role: Deep muscular holding and pelvic or hip reconstructive exposure · Price (USD): \$44.00 / \$49.00
1	Protective Sleeve / Box	Available Lengths & Handles: Standard Transit Packaging · Primary Clinical Role: — · Price (USD): Included

PRODUCT DESCRIPTION

The Kocher Langenbeck Retractor is an indispensable handheld surgical retractor recognized worldwide for providing safe, non-traumatic soft-tissue and muscular retraction during general, orthopedic, and plastic surgery. Designed with an elongated right-angled blade terminating in a downward retaining lip, this instrument allows surgeons to retract muscle planes and wound margins securely without slipping or causing focal tissue ischemia.

Built for complete procedural adaptability, this Kocher Langenbeck Retractor series is handcrafted from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) across two overall instrument lengths: 20.3 cm (\$44) and 22 cm (\$49), with a 215 mm standard frame text profile. Surgeons can customize the instrument across three shaft lengths (10 cm, 12 cm, and 14 cm) paired with either a classic Finger Ring Handle or a lightweight Hollow Grip Handle engineered to enhance comfort and reduce fatigue during long surgeries. The extensive blade layout encompasses 20 standardized dimensions from 15 × 85 mm to 80 × 16 mm.

Forged with an anti-glare satin matte finish, Retractor Maker manufactures the Kocher Langenbeck Retractor to withstand repetitive hospital autoclave cycles, providing operating theaters worldwide with exceptional durability, corrosion resistance, and complete clinical dependability.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)