



RETRACTOR MAKER

PRODUCT BROCHURE OF

W.M. Lemmon Deaver Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



W.M. Lemmon Deaver Retractor

SKU RM-WMLDR-395-1

\$119 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
25 mm Wide Blade — 25 mm wide x 235 mm working surface	RM-WMLDR-395-1-1	\$119
50 mm Wide Blade — 50 mm wide x 235 mm working surface	RM-WMLDR-395-1-2	\$119

PRODUCT DESCRIPTION

The W.M. Lemmon Deaver Retractor is a premier deep-cavity surgical retractor engineered specifically to provide safe, broad-surface tissue displacement of abdominal viscera, lung lobes, and deep pelvic organs during major thoracic, laparotomy, and pelvic procedures. Designed with an elongated, sweeping curved blade, this instrument allows surgical assistants to maintain constant counter-traction across heavy visceral structures without causing localized parenchymal contusions or tears.

Built to withstand high-frequency autoclave sterilization cycles without degrading, this W.M. Lemmon Deaver Retractor configuration features an extraordinary overall tool length of 39.5 cm (\$119) that provides reaching clearance into the deepest anatomical corridors. Handcrafted from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic), the active blade length measures 235 mm along its working surface. The instrument offers two selectable blade width options: a 25 mm wide blade layout for focused visceral holding and deep vascular exposure, and an expansive 50 mm wide blade layout for broad liver, splenic, or intestinal displacement. Its ergonomic flat handle delivers fatigue-free manual leverage.

Forged with an anti-glare satin matte finish, Retractor Maker produces the W.M. Lemmon Deaver Retractor to uphold strict clinical tolerances, delivering exceptional structural rigidity, corrosion resistance, and complete autoclave durability for general surgery operating theaters, trauma centers, and cardiothoracic suites worldwide.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)