



RETRACTOR MAKER

PRODUCT BROCHURE OF

Bookwalter Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Bookwalter Retractor

SKU RM-BR-1

\$1,999 USD

COMPONENT LIST

QUANTITY	COMPONENT	DIMENSIONS / DESCRIPTION
1	Bookwalter Field Post	12" (30.5 cm)
1	Bookwalter Horizontal Bar	Standard
1	Bookwalter Post Coupler	Standard
1	Bookwalter Horizontal Flex Bar	1" Diameter (2.5 cm)
2	Bookwalter Standard Ratchets	Non-Tilt
5	Bookwalter Tilt Ratchet Mechanisms	Multi-Angle
1	Bookwalter Oval Ring, Medium	10½" × 14½" (26.7 × 36.8 cm)
1	Bookwalter Segmented Ring, Medium	10½" × 19" (26.7 × 48.3 cm)
2	Bookwalter Balfour Blades	2⅝" × 3" (6.7 × 7.6 cm)
1	Bookwalter Kelly Blade	2" × 3" (5.1 × 7.6 cm)
1	Bookwalter Kelly Blade	2" × 4" (5.1 × 10.2 cm)
1	Bookwalter Kelly Blade	2" × 5" (5.1 × 12.7 cm)
1	Bookwalter Kelly Blade	2" × 6" (5.1 × 15.2 cm)
1	Bookwalter Malleable Blade	3" × 6" (7.6 × 15.2 cm)

QUANTITY	COMPONENT	DIMENSIONS / DESCRIPTION
2	Bookwalter Malleable Blades	2" × 6" (5.1 × 15.2 cm)
1	Bookwalter Harrington Blade	1" × 5 ⁵ / ₈ " × 2 ¹ / ₂ " (2.5 × 14.3 × 6.4 cm)
1	Wooden Storage Box	Fitted Protective Case

PRODUCT DESCRIPTION

The Bookwalter Retractor is a table-mounted, self-retaining retractor system. Instead of relying on assistants to hold individual retractors, the set anchors to the operating table and holds each blade in a fixed position, delivering steady, hands-free exposure and surgical field visibility during major open abdominal procedures such as nephrectomy, prostatectomy and laparotomy.

The mounting framework consists of a 12" (30.5 cm) Bookwalter field post, a horizontal bar, a post coupler and a 1" (2.5 cm) diameter horizontal flex bar. This assembly supports a medium oval ring measuring 10¹/₂" × 14¹/₂" (26.7 × 36.8 cm) and a medium segmented ring measuring 10¹/₂" × 19" (26.7 × 48.3 cm), allowing the retracting frame to be positioned over the incision and adjusted to the working area.

Two standard non-tilt ratchets and five multi-angle tilt ratchet mechanisms clamp to the ring and hold each blade at the selected depth and angle. The included blades comprise two Balfour blades, four Kelly blades in 3", 4", 5" and 6" lengths, three malleable blades and a Harrington blade, so exposure can be adapted to the depth and shape of the incision.

All components are manufactured from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic). The Bookwalter Retractor is brand new, reusable and fully autoclavable, and is supplied non-sterile — it must be cleaned and sterilized according to the healthcare facility's validated procedures before use. Components are securely arranged in a dedicated, perforated wooden box for organized storage and processing. Retractor Maker manufactures the system in Sialkot, Pakistan and supplies it directly to hospitals, surgeons and distributors, with a minimum order quantity of one full set. The Bookwalter Retractor is manufactured, inspected and supplied directly by Retractor Maker.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)