



RETRACTOR MAKER

PRODUCT BROCHURE OF

Tubular Retractor System

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved




SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Tubular Retractor System

SKU RM-TRS-MIS-1

\$2,199 USD

COMPONENT LIST

QUANTITY	COMPONENT	DIMENSIONS / DESCRIPTION
1	MIS Tubular Retractor Neupsy Key	Diameter 12-14 mm, length 30 mm
1	MIS Tubular Retractor Neupsy Key	Diameter 12-14 mm, length 40 mm
1	MIS Tubular Retractor Neupsy Key	Diameter 12-14 mm, length 50 mm
1	MIS Tubular Retractor Neupsy Key	Diameter 12-14 mm, length 60 mm
1	MIS Tubular Retractor Neupsy Key	Diameter 12-14 mm, length 70 mm
1	MIS Tubular Retractor Neupsy Key	Diameter 12-14 mm, length 80 mm
1	MIS Tubular Retractor Neupsy Key	Diameter 12-14 mm, length 90 mm
1	MIS Tubular Retractor Neupsy Key	Diameter 14-16 mm, length 30 mm
1	MIS Tubular Retractor Neupsy Key	Diameter 14-16 mm, length 40 mm
1	MIS Tubular Retractor Neupsy Key	Diameter 14-16 mm, length 50 mm
1	MIS Tubular Retractor Neupsy Key	Diameter 14-16 mm, length 60 mm
1	MIS Tubular Retractor Neupsy Key	Diameter 14-16 mm, length 70 mm
1	MIS Tubular Retractor Neupsy Key	Diameter 14-16 mm, length 80 mm
1	MIS Tubular Retractor Neupsy Key	Diameter 14-16 mm, length 90 mm

QUANTITY	COMPONENT	DIMENSIONS / DESCRIPTION
1	MIS Tubular Retractor Neupsy Key	Diameter 16-18 mm, length 30 mm
1	MIS Tubular Retractor Neupsy Key	Diameter 16-18 mm, length 40 mm
1	MIS Tubular Retractor Neupsy Key	Diameter 16-18 mm, length 50 mm
1	MIS Tubular Retractor Neupsy Key	Diameter 16-18 mm, length 60 mm
1	MIS Tubular Retractor Neupsy Key	Diameter 16-18 mm, length 70 mm
1	MIS Tubular Retractor Neupsy Key	Diameter 16-18 mm, length 80 mm
1	MIS Tubular Retractor Neupsy Key	Diameter 16-18 mm, length 90 mm
1	MIS Tubular Retractor Neupsy Key	Diameter 18-20 mm, length 30 mm
1	MIS Tubular Retractor Neupsy Key	Diameter 18-20 mm, length 40 mm
1	MIS Tubular Retractor Neupsy Key	Diameter 18-20 mm, length 50 mm
1	MIS Tubular Retractor Neupsy Key	Diameter 18-20 mm, length 60 mm
1	MIS Tubular Retractor Neupsy Key	Diameter 18-20 mm, length 70 mm
1	MIS Tubular Retractor Neupsy Key	Diameter 18-20 mm, length 80 mm
1	MIS Tubular Retractor Neupsy Key	Diameter 18-20 mm, length 90 mm
1	MIS Tubular Retractor Neupsy Key	Diameter 27-29 mm, length 30 mm
1	MIS Tubular Retractor Neupsy Key	Diameter 27-29 mm, length 40 mm
1	MIS Tubular Retractor Neupsy Key	Diameter 27-29 mm, length 50 mm
1	MIS Tubular Retractor Neupsy Key	Diameter 27-29 mm, length 60 mm
1	MIS Tubular Retractor Neupsy Key	Diameter 27-29 mm, length 70 mm
1	MIS Tubular Retractor Neupsy Key	Diameter 27-29 mm, length 80 mm
1	MIS Tubular Retractor Neupsy Key	Diameter 27-29 mm, length 90 mm
1	MIS Tubular Retractor Dilator	Diameter 6 mm
1	MIS Tubular Retractor Dilator	Diameter 8 mm
1	MIS Tubular Retractor Dilator	Diameter 10 mm
1	MIS Tubular Retractor Dilator	Diameter 12 mm
1	MIS Tubular Retractor Dilator	Diameter 14 mm

QUANTITY	COMPONENT	DIMENSIONS / DESCRIPTION
1	MIS Tubular Retractor Dilator	Diameter 16 mm
1	MIS Tubular Retractor Dilator	Diameter 18 mm
1	Tubular Flexible Arm	Multi-jointed articulating holding arm
1	Table Clamp Set	Operating table rail mount
1	Sterilization Box	Fitted case for organized reprocessing

PRODUCT DESCRIPTION

The Tubular Retractor System is an essential minimally invasive surgical instrument platform engineered to establish muscle-sparing access corridors during complex spinal interventions. By facilitating sequential tissue dilation rather than traumatic muscle stripping, this system dramatically minimizes tissue disruption, reduces postoperative pain, and shortens patient recovery timelines.

Engineered for procedural versatility across both cervical and lumbar spine segments, the Tubular Retractor System features an extensive complement of specialized tubular ports spanning 12-14 mm, 14-16 mm, 16-18 mm, 18-20 mm, and 27-29 mm diameters across lengths from 30 mm to 90 mm. The set includes sequential dilators from 6 mm to 18 mm, a rigid table rail clamp, and a flexible articulating arm that anchors securely to the operating table to lock the working tube rigidly over the target disc space.

Precision-machined from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with a reflection-free satin finish, Retractor Maker supplies the complete Tubular Retractor System non-sterile and fully autoclavable inside a dedicated sterilization box, ensuring optimal sterility management and organized OR turnover.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)