



RETRACTOR MAKER

PRODUCT BROCHURE OF

Jansen Mastoid Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved




SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Jansen Mastoid Retractor

SKU RM-JMR-1

\$69 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
3"	RM-JMR-1-1	\$69
3 ¼"	RM-JMR-1-2	\$74
3 ⅝"	RM-JMR-1-3	\$79
4"	RM-JMR-1-4	\$84
4 ¼"	RM-JMR-1-5	\$89
4 ½"	RM-JMR-1-6	\$94

PRODUCT DESCRIPTION

The Jansen Mastoid Retractor is an internationally acclaimed self-retaining surgical retractor engineered specifically to provide stable, hands-free soft-tissue and periosteal retraction in confined otologic, scalp, and microsurgical fields. Designed with outward-curving multi-prong blades, this instrument holds postauricular wound margins open under calibrated mechanical tension, eliminating assistant hand fatigue during delicate mastoid drilling.

Built for complete procedural adaptability, this Jansen Mastoid Retractor series is highly customizable across six overall tool lengths: 3" (\$69), 3 ¼" (\$74), 3 ⅝" (\$79), 4" (\$84), 4 ¼" (\$89), and 4 ½" (\$94). Surgeons can select between 3×3 or 4×4 prong arrays in Sharp tips for firm fascial purchase or Blunt tips to prevent local tissue trauma, all standardized at an 8 mm functional prong depth. The instrument offers spread configurations including 15 mm blade length × 70 mm spread or 20 mm blade length × 64 mm spread, driven by a fine-pitch thumbscrew locking mechanism.

Machined from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with a reflection-free satin matte finish, Retractor Maker produces the Jansen Mastoid Retractor to strict clinical tolerances, delivering high corrosion resistance, smooth mechanical locking, and complete autoclave durability for ENT and neurosurgical suites worldwide.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)