



RETRACTOR MAKER

PRODUCT BROCHURE OF

Brewster Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Brewster Retractor

SKU RM-BRW-1

\$159 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
Class 'E' (1.9 × 3.8 cm)	RM-BRW-1-1	\$159
Class 'D' (4.4 × 2.5 cm)	RM-BRW-1-2	\$164
Class 'C' (5.1 × 3.5 cm)	RM-BRW-1-3	\$169
Class 'B' (5.7 × 4.4 cm)	RM-BRW-1-4	\$174
Class 'A' (6.3 × 4.7 cm)	RM-BRW-1-5	\$179

PRODUCT DESCRIPTION

The Brewster Retractor is an adaptable handheld surgical retractor engineered specifically to provide reliable, non-slip exposure of deep soft-tissue layers, muscle planes, and wound borders. Designed with a robust shank and an anatomically contoured distal blade, this instrument allows surgical assistants to apply continuous, firm counter-traction without crushing fragile tissue margins.

Built for surgical versatility across diverse anatomical depths, the Brewster Retractor series features five distinct standardized blade size classes: Size 'A' (6.3 × 4.7 cm, \$179), Size 'B' (5.7 × 4.4 cm, \$174), Size 'C' (5.1 × 3.5 cm, \$169), Size 'D' (4.4 × 2.5 cm, \$164), and Size 'E' (1.9 × 3.8 cm, \$159). Each size integrates a smooth retaining lip that hooks securely behind fascia and muscle margins, preventing blade slippage while maintaining a clear operative view.

Manufactured from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with a reflection-free satin matte finish, Retractor Maker produces the Brewster Retractor to withstand heavy mechanical traction and repetitive hospital sterilization cycles, providing general, plastic, and orthopedic surgical teams with long-lasting

clinical dependability.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)