



RETRACTOR MAKER

PRODUCT BROCHURE OF

Fomon Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Fomon Retractor

SKU RM-FMR-17-1

\$44 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
Standard Configuration — 17 cm length, 5 × 10 mm ball tip, stainless steel	RM-FMR-17-1-1	\$44
Insulated Layout — 17 cm length, 5 × 10 mm ball tip, dielectric coating	RM-FMR-17-1-2	\$54

PRODUCT DESCRIPTION

The Fomon Retractor is an internationally acclaimed rhinoplasty and facial plastic surgical instrument engineered specifically to provide atraumatic soft-tissue retraction and clear visual access within compact anatomical corridors. Designed with an ultra-smooth ball-shaped tip, this retractor hooks gently beneath the alar cartilage and mucosal linings, shielding delicate tissues from sharp puncturing or focal pressure trauma.

Built for optimal tactile control during nasal remodeling, this Fomon Retractor model features an overall tool length of 17 cm (6.7 inches). Its working end incorporates a dedicated 5 mm × 10 mm ball-tip measurement contoured to distribute counter-traction evenly across the lower lateral cartilage and alar rim. Surgeons can choose between a Standard Stainless Steel layout (\$44) and an Insulated layout (\$54) coated with a durable non-conductive polymer to protect adjacent skin and mucosal margins during active electrosurgery.

Forged from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with a reflection-free satin matte finish, Retractor Maker manufactures the Fomon Retractor to strict clinical standards, delivering high corrosion resistance, unyielding structural rigidity, and complete autoclave durability for cosmetic, ENT, and reconstructive clinics worldwide.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)