



RETRACTOR MAKER

PRODUCT BROCHURE OF

Lateral Fat Pad Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Lateral Fat Pad Retractor

SKU RM-LFPR-10IN-1

\$84 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
Lateral Fat Pad Retractor	Dimensions / Details: 10" (25.4 cm) overall length, 41 mm wide polished curved blade · Quantity Included: 1	\$84
Protective Sleeve / Box	Dimensions / Details: Individual transit packaging · Quantity Included: 1	\$84

PRODUCT DESCRIPTION

The Lateral Fat Pad Retractor is an indispensable orthopedic knee surgical instrument engineered specifically to provide safe, non-traumatic displacement of the infrapatellar fat pad and surrounding synovial tissues during total knee arthroplasty. Designed with a wide, curved blade, this retractor allows surgical teams to hold adipose structures cleanly away from the lateral compartment, exposing the tibial plateau and lateral meniscus for resection.

Built to maximize tactile feedback and procedural efficiency, this Lateral Fat Pad Retractor configuration features an overall tool length of 10 inches (25.4 cm, \$84). Handcrafted from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic), its working end incorporates a 41 mm wide curved blade profile with a smooth, polished working surface specifically engineered to minimize tissue friction, preventing shredding of delicate adipose lobules. The heavy-duty handle provides operating assistants with steady manual leverage throughout prolonged joint preparation.

Forged with an anti-glare satin matte finish, Retractor Maker produces the Lateral Fat Pad Retractor to withstand heavy mechanical loads and repeated hospital autoclave cycles, providing knee arthroplasty and orthopedic reconstructive teams worldwide with exceptional durability, corrosion resistance, and complete clinical dependability.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)