



RETRACTOR MAKER

PRODUCT BROCHURE OF

Hoer Nerve Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved




SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Hoen Nerve Retractor

SKU RM-HNR-22-1

\$24 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
Hoen Nerve Retractor	22 cm overall length, curved atraumatic nerve tip	\$24
Protective Sleeve / Box	Individual transit packaging	\$24

PRODUCT DESCRIPTION

The Hoen Nerve Retractor is an internationally acclaimed neurosurgical and spinal instrument engineered specifically to provide safe, non-traumatic nerve root and thecal sac displacement during spinal decompression procedures. Designed with an ultra-smooth, curved distal tip, this retractor allows surgeons to gently explore, mobilize, and shield delicate neural fibers without causing localized traction neuropraxia or dural lacerations.

Built for delicate nerve root displacement in deep operative corridors, this Hoen Nerve Retractor configuration features an overall tool length of 22 cm (\$24). Handcrafted from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) to offer absolute corrosion resistance and long life, its slender, balanced shaft provides the operating surgeon with refined tactile sensitivity. The knurled ergonomic handle ensures non-slip fingertip control under high-magnification operating microscopes.

Manufactured to strict clinical tolerances with an anti-glare satin matte finish, Retractor Maker supplies the Hoen Nerve Retractor non-sterile and fully autoclavable inside protective packaging, providing neurosurgeons and spine specialists worldwide with unmatched micro-dissection control and clinical dependability.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)