



RETRACTOR MAKER

PRODUCT BROCHURE OF

Clodius Orbital Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Clodius Orbital Retractor

SKU RM-COR-1

\$24 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
Clodius Orbital Retractor	Handheld periorbital and blepharoplasty retractor	\$24
Protective Sleeve / Box	Individual transit packaging	\$24

PRODUCT DESCRIPTION

The Clodius Orbital Retractor is a specialized periorbital surgical retractor engineered specifically to provide safe, non-traumatic retraction of orbital fat pads and delicate cutaneous flaps during aesthetic blepharoplasty and orbital reconstruction. Designed with an ultra-smooth, contoured blade, this instrument shields the ocular globe while establishing a clean visual field for fat excision and septum repositioning.

Built for precision microsurgical clearance, this Clodius Orbital Retractor configuration (\$24) features a slender, balanced shaft that provides the operating surgeon with refined tactile sensitivity in confined ocular spaces. Its curved spatula tip incorporates rounded, non-traumatic margins that seat securely against the orbital rim, preventing accidental globe compression or vascular bruising. The ergonomic flat handle ensures fatigue-free control during delicate eyelid remodeling.

Forged from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with a reflection-free satin matte finish, Retractor Maker manufactures the Clodius Orbital Retractor to withstand repetitive hospital autoclave cycles, providing oculoplastic, cosmetic, and ophthalmic surgical suites worldwide with high corrosion resistance and clinical precision.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)