



RETRACTOR MAKER

PRODUCT BROCHURE OF

Holscher Retractors

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Holscher Retractors

SKU RM-HLSR-1

\$43 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
7 mm Wide — 9 ½" (24.13 cm)	RM-HLSR-1-1	\$43
9 mm Wide — 9 ½" (24.13 cm)	RM-HLSR-1-2	\$46
13 mm Wide — 9 ½" (24.13 cm)	RM-HLSR-1-3	\$49

PRODUCT DESCRIPTION

Holscher Retractors are specialized spinal microsurgery instruments engineered specifically to provide safe, non-traumatic retraction and protection of the thecal sac and spinal nerve roots during posterior lumbar decompression procedures. Designed with an elongated slender shaft, these retractors allow spine surgeons to work comfortably in deep operative corridors while keeping manual handles safely outside the microscopic line of sight.

Built to fit varying anatomical requirements, the Holscher Retractors series maintains a uniform overall tool length of 9 ½ inches (24.13 cm) that provides exceptional reaching depth during disc excisions. The line is manufactured across three selectable working end widths: 7 mm (\$43) for narrow interlaminar windows, 9 mm (\$46) for standard nerve root mobilization, and 13 mm (\$49) for broad thecal sac protection. Its smoothly contoured blade edges cradle delicate nerve fibers securely without causing localized traction neuropraxia.

Forged from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with a reflection-free satin matte finish, Retractor Maker produces Holscher Retractors to withstand repetitive hospital autoclave cycles, providing spinal, neurosurgical, and orthopedic operating theaters worldwide with unyielding structural rigidity, high corrosion resistance, and complete clinical durability.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)