



RETRACTOR MAKER

PRODUCT BROCHURE OF

Wullstein Weitlaner Self Retaining Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Wullstein Weitlaner Self Retaining Retractor

SKU RM-WWSRR-1

\$44 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
Compact (4", 4 ¼", 4 ½") — Straight / Curved; 2×2 or 2×3 prongs	RM-WWSRR-1-1	\$44
Intermediate (5 ⅞", 5.5", 6", 6 ½") — Straight / Curved; 2×3, 3×3, or 3×4 prongs	RM-WWSRR-1-2	\$59
Extended (8", 9 ½", 10") — Straight / Curved; 3×3 or 3×4 prongs	RM-WWSRR-1-3	\$79

PRODUCT DESCRIPTION

The Wullstein Weitlaner Self Retaining Retractor is an internationally acclaimed self-retaining surgical spreader engineered specifically to provide dependable, hands-free soft-tissue and periosteal retraction in confined surgical corridors. Designed with outward-curving multi-prong jaws and an integrated mechanical ratchet lock, this instrument maintains constant, calibrated counter-traction across opposing wound margins, eliminating assistant hand fatigue while preserving an open visual field under operating microscopes.

Built to withstand frequent sterilization and autoclaving protocols, this Wullstein Weitlaner Self Retaining Retractor series features an extensive sizing matrix across ten overall lengths: 4" (\$44), 4 ¼" (\$49), 4 ½" (\$54), 5 ⅞" (\$59), 5.5" (\$64), 6" (\$69), 6 ½" (\$74), 8" (\$79), 9 ½" (\$84), and 10" (\$89). Handcrafted from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic), the instrument offers Straight or Curved jaw arm profiles paired with 2×2, 2×3, 3×3, or 3×4 prong density options. Surgeons can select Blunt tips optimized for safe, extra-traumatic retraction near sensitive vascular pathways, or Sharp tips engineered to deliver an unyielding, non-slip purchase on tougher dermal or fascial layers.

Forged with an anti-glare satin matte finish, Retractor Maker manufactures the Wullstein Weitlaner Self Retaining Retractor to uphold strict clinical standards, providing ENT suites, neurosurgical operating theaters, and plastic reconstructive clinics worldwide with exceptional corrosion resistance, mechanical stability, and complete autoclave durability.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)