



RETRACTOR MAKER

PRODUCT BROCHURE OF

Williams Distal Radius Fracture Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Williams Distal Radius Fracture Retractor

SKU RM-WDRFR-45-1

\$189 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
Curved Left Blade — 4.5" (11.43 cm) length, blunt angled prongs, left curve	RM-WDRFR-45-1-1	\$189
Curved Right Blade — 4.5" (11.43 cm) length, blunt angled prongs, right curve	RM-WDRFR-45-1-2	\$189

PRODUCT DESCRIPTION

The Williams Distal Radius Fracture Retractor is an indispensable orthopedic hand surgical instrument engineered specifically to provide safe, slip-resistant bone leverage and soft-tissue displacement during volar plating of distal radius fractures. Designed to wrap around the radial cortex through the Henry or modified flexor carpi radialis (FCR) approach, this retractor holds the pronator quadratus and flexor tendons aside cleanly without compressing the median nerve.

Built to withstand high-frequency autoclave sterilization cycles without degrading, this Williams Distal Radius Fracture Retractor configuration features an overall tool length of 4.5 inches (11.43 cm, \$189). Handcrafted from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic), the prong profile configuration features blunt left and right angled prongs specifically contoured to maintain secure, slip-resistant bone contact on the volar radius border. Surgeons can choose between two distinct geometric profiles: a Curved Left blade layout or a Curved Right blade layout to match anatomical laterality. Its ergonomic handle delivers fatigue-free manual leverage.

Forged with an anti-glare satin matte finish, Retractor Maker manufactures the Williams Distal Radius Fracture Retractor to uphold strict clinical tolerances, providing hand surgery centers, trauma units, and orthopedic reconstructive suites worldwide with exceptional durability, corrosion resistance, and complete autoclave durability.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractor maker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractor maker.com
retractor manufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)