



RETRACTOR MAKER

PRODUCT BROCHURE OF

Kristeller Vaginal Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Kristeller Vaginal Retractor

SKU RM-KVR-1

\$44 USD

COMPONENT LIST

QUANTITY	COMPONENT	DIMENSIONS / DESCRIPTION
1	Narrow Blades (75×10, 110×14, 110×17 mm)	Available Lengths: 18 cm or 22 cm · Primary Clinical Role: Pediatric, nulliparous, and narrow vaginal canal exposure · Price (USD): \$44.00 / \$49.00
1	Medium Blades (75×23, 80×25, 115×23 mm)	Available Lengths: 18 cm or 22 cm · Primary Clinical Role: Standard adult vaginal examination and minor pelvic surgery · Price (USD): \$44.00 / \$49.00
1	Broad/Deep Blades (100×32, 115×26, 115×29, 115×32 mm)	Available Lengths: 18 cm or 22 cm · Primary Clinical Role: Deep pelvic floor repair and vaginal hysterectomy in parous patients · Price (USD): \$44.00 / \$49.00
1	Protective Sleeve / Box	Available Lengths: Standard Transit Packaging · Primary Clinical Role: — · Price (USD): Included

PRODUCT DESCRIPTION

The Kristeller Vaginal Retractor is a foundational gynecological surgical instrument engineered specifically to provide broad, atraumatic displacement of the vaginal walls and pelvic floor during diagnostic and operative pelvic interventions. Designed with a smoothly contoured scoop blade, this instrument depresses vaginal mucosa cleanly, establishing an expansive visual corridor to the cervix and fornices without causing mucosal laceration.

Built for complete clinical adaptability, this Kristeller Vaginal Retractor series is handcrafted from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) across two overall tool lengths: 18 cm (\$44) and 22 cm (\$49). Surgeons can choose from an extensive matrix of ten blade footprint selections spanning lengths from 75

mm to 115 mm and widths from 10 mm to 32 mm: 75 × 10 mm, 75 × 23 mm, 80 × 25 mm, 100 × 32 mm, 110 × 14 mm, 110 × 17 mm, 115 × 23 mm, 115 × 26 mm, 115 × 29 mm, and 115 × 32 mm. Its ergonomic drop handle delivers steady manual leverage for operating assistants.

Forged with an anti-glare satin matte finish, Retractor Maker manufactures the Kristeller Vaginal Retractor to withstand repetitive hospital sterilization cycles, providing obstetrical, gynecological, and urogynecological teams worldwide with exceptional durability, corrosion resistance, and complete autoclave durability.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)