



RETRACTOR MAKER

PRODUCT BROCHURE OF

Iron Intern Single Arm

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Iron Intern Single Arm

SKU RM-IISA-SERIES-1

\$999 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
17" Standard Model	RM-IISA-1-1	\$999
22" Standard Model	RM-IISA-2-1	\$1,199
22" Jumbo Heavy Duty	RM-IISA-3-1	\$1,299
22" Duplex Variant	RM-IISA-4-1	\$1,399
22" Collapsible Variant	RM-IISA-5-1	\$1,499
26" Standard Model	RM-IISA-6-1	\$1,599
26" Jumbo Heavy Duty	RM-IISA-7-1	\$1,699

PRODUCT DESCRIPTION

The Iron Intern Single Arm is an internationally acclaimed table-mounted surgical retractor system engineered specifically to provide rock-solid, hands-free visceral retraction and structural tissue displacement during open general, bariatric, and laparoscopic procedures. Designed to clamp securely onto standard operating table side rails, this mechanical holding arm replaces manual assistant fatigue with unwavering static counter-traction, maintaining a wide, stationary operative corridor throughout lengthy abdominal and pelvic operations.

Built to clinical framework specifications, this Iron Intern Single Arm series features an extensive selection matrix spanning three longitudinal pole reaches and specialized duty tiers. Handcrafted from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic), configurations include a 17-inch (43 cm) Standard arm (\$999) for pediatric and compact fields; a 22-inch (56 cm) Standard arm (\$1,199); a 22-inch Jumbo Heavy Duty arm

(\$1,299) for bariatric loads; a 22-inch Duplex Variant (\$1,399) supporting multiple accessory clamps; a 22-inch Collapsible Variant (\$1,499) for quick clearance without releasing table clamps; a 26-inch (66 cm) Standard arm (\$1,599) for deep open laparotomies; and an ultra-rigid 26-inch Jumbo Heavy Duty arm (\$1,699) built for complex liver transplants and extensive oncology.

Forged with an anti-glare satin matte finish, Retractor Maker manufactures the Iron Intern Single Arm to withstand repetitive hospital autoclave cycles, providing general surgery suites, bariatric centers, and trauma operating theaters worldwide with exceptional structural rigidity, corrosion resistance, and complete autoclave durability.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractor maker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractor maker.com
retractor manufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)