



RETRACTOR MAKER

PRODUCT BROCHURE OF

Richardson Eastman Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Richardson Eastman Retractor

SKU RM-RER-SERIES-1

\$54 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
10 ¼" (approx. 26.0 cm) — Dual right-angled concave blades	RM-RER-SERIES-1-1	\$54
10 ½" (approx. 26.7 cm) — Dual right-angled concave blades	RM-RER-SERIES-1-2	\$59
10 ¾" (approx. 27.3 cm) — Dual right-angled concave blades	RM-RER-SERIES-1-3	\$64
11 ¼" (approx. 28.6 cm) — Dual right-angled concave blades	RM-RER-SERIES-1-4	\$69

PRODUCT DESCRIPTION

The Richardson Eastman Retractor is an internationally acclaimed double-ended surgical retractor engineered specifically to provide dependable, multi-depth wound margin displacement and abdominal wall exposure during major laparotomies and pelvic procedures. Designed with two distinct right-angled concave blades on opposite ends of a contoured central handle, this instrument allows surgical teams to transition effortlessly between shallow subcutaneous tissue retraction and deep muscular holding without switching instruments.

Built to provide exceptional durability and corrosion resistance, this Richardson Eastman Retractor series is available across four overall lengths: 10 ¼ inches (\$54), 10 ½ inches (\$59), 10 ¾ inches (\$64), and 11 ¼ inches (\$69). Handcrafted from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic), the blade geometry layout incorporates an asymmetric double-ended design combining two distinct right-angled concave blades with smooth, rounded retaining lips that prevent tissue slippage or localized pressure contusions. Available to order as a Single instrument or as a Set of Two retractors, its central handle ensures fatigue-free manual leverage.

Forged with an anti-glare satin matte finish, Retractor Maker manufactures the Richardson Eastman Retractor to withstand repetitive hospital autoclave cycles, delivering high corrosion resistance, unyielding structural rigidity, and complete autoclave durability for general surgery suites, obstetrical theaters, and trauma centers worldwide.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractor maker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractor maker.com
retractor manufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)