



RETRACTOR MAKER

PRODUCT BROCHURE OF

Finochietto-Baby Rib Spreader

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.

Finochietto-Baby Rib Spreader

SKU RM-FBRS-100-1

\$299 USD



COMPONENT LIST

QUANTITY	COMPONENT	DIMENSIONS / DESCRIPTION
1	Parallel Slide Guide Rail Track Frame	Miniature track system designed to slide open within tight pediatric spaces
2 (Fixed Pair)	Fixed Angled Blades	31 mm × 31 mm bone-gripping lateral blocks with circular fenestrations
1	Integrated Turning Screw Lock	Manual structural lock positioned to securely hold frame expansion
2	Ergonomic Forceps Extension Arms	Ultra-lightweight structured handles engineered with smooth edges
1	Protective Storage / Box	Standard transit packaging

PRODUCT DESCRIPTION

The Finochietto-Baby Rib Spreader is a specialized pediatric thoracic surgical retractor engineered specifically to provide safe, gradual rib separation and chest wall exposure during infant thoracotomies and pediatric cardiothoracic interventions. Designed with a scaled-down slide guide rail, this instrument allows surgeons to distribute mechanical spreading forces gently across fragile ribs, preventing costal fractures while maintaining an unobstructed visual corridor into the thoracic cavity.

Built specifically for tight pediatric operative fields, this Finochietto-Baby Rib Spreader configuration (\$299) features fixed fenestrated lateral blades measuring 31 mm × 31 mm that are anatomically optimized for neonatal and pediatric anatomy. Handcrafted from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic), the parallel slide track opens to a 100 mm maximum expansion capability. The circular blade fenestrations ensure positive purchase on delicate rib borders, while the integrated turning screw lock securely retains the

chosen retraction spread without slipping, eliminating assistant fatigue during delicate procedures.

Forged with an anti-glare satin matte finish, Retractor Maker produces the Finochietto-Baby Rib Spreader to withstand repetitive hospital sterilization cycles, providing pediatric cardiac operating rooms, children's hospitals, and congenital thoracic centers worldwide with exceptional structural durability, corrosion resistance, and complete autoclave dependability.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractor maker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractor maker.com
retractor manufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)