



RETRACTOR MAKER

PRODUCT BROCHURE OF

Curtin Skin and Tissue Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Curtin Skin and Tissue Retractor

SKU RM-CSTR-SKIN-1

\$49 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
Curtin Skin and Tissue Retractor	14 cm overall length, 15.0 mm wide head, 6 sharp prongs	\$49
Protective Sleeve / Box	Individual transit packaging	\$49

PRODUCT DESCRIPTION

The Curtin Skin and Tissue Retractor is a specialized handheld dermatologic and plastic surgery retractor engineered specifically to provide firm, non-slip skin flap elevation and subcutaneous tissue exposure. Designed with an expansive multi-prong rake head, this instrument allows surgeons to apply uniform manual counter-traction across broad dermal layers without tearing fragile cutaneous margins.

Built for steady manual leverage in superficial operative fields, this Curtin Skin and Tissue Retractor configuration features an overall length of 14 cm (\$49) and a 15.0 mm wide working end. Its working end is formatted with 6 sharp prongs structured to provide precise tissue engagement, distributing mechanical tension evenly to prevent tissue slippage during sharp dissection. The low-profile ergonomic handle ensures a stable, fatigue-free grip for operating assistants during lengthy aesthetic procedures.

Manufactured from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with a reflection-free satin matte finish, Retractor Maker produces the Curtin Skin and Tissue Retractor to withstand repetitive hospital sterilization cycles, providing plastic, reconstructive, and dermatologic surgical teams with high corrosion resistance, unyielding precision, and complete autoclave durability.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)