



RETRACTOR MAKER

PRODUCT BROCHURE OF

Knapp Self-Retaining Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Knapp Self-Retaining Retractor

SKU RM-KSRR-10-1

\$79 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
Knapp Self-Retaining Retractor	Dimensions / Details: 10 cm overall length, 26 mm width, 4×4 prongs (Straight/Angled, Sharp/Blunt) · Primary Clinical Role: Self-retaining lacrimal sac and micro-incision exposure · Price (USD): \$79.00	\$79
Protective Sleeve / Box	Dimensions / Details: Standard Transit Packaging · Primary Clinical Role: — · Price (USD): Included	\$79

PRODUCT DESCRIPTION

The Knapp Self-Retaining Retractor is a specialized ophthalmic and micro-surgical spreader engineered specifically to provide reliable, hands-free soft-tissue and wound margin exposure during delicate facial, ocular, and hand procedures. Designed with outward-curving prongs, this instrument holds small incisions open under continuous mechanical tension, eliminating assistant hand fatigue while maintaining a static operative corridor.

Built to clinical framework standards, this Knapp Self-Retaining Retractor maintains a uniform overall instrument length of 10 cm (\$79) paired with a 26 mm working blade width. Handcrafted from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic), the instrument incorporates a 4×4 outward-curving prong array selectable in either a Straight configuration layout or an Angled configuration layout to bypass orbital margins. Surgeons can choose between Blunt tips designed for a gentle yet effective grip on tissues or Sharp tips designed for highly secure tissue engagement in dense periosteum.

Forged with an anti-glare satin matte finish, Retractor Maker manufactures the Knapp Self-Retaining Retractor to withstand repetitive hospital autoclave cycles, providing plastic, ophthalmic, and ENT surgeons worldwide

with exceptional mechanical stability, corrosion resistance, and complete clinical durability.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)