



RETRACTOR MAKER

PRODUCT BROCHURE OF

Bishop Lip And Cheek Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.

Donald Cameron



Bishop Lip And Cheek Retractor

SKU RM-BLCR-1

\$24 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
Bishop Lip And Cheek Retractor	13.5 cm overall length, 27 mm double-ended blades	\$24
Protective Sleeve / Box	Individual transit packaging	\$24

PRODUCT DESCRIPTION

The Bishop Lip And Cheek Retractor is a versatile, double-ended surgical retractor engineered specifically to provide atraumatic displacement of the cheeks, lips, and oral mucosa during dental, maxillofacial, and cleft palate procedures. Designed with smooth, rounded margins, this instrument holds intraoral soft tissues away from surgical drills and suture lines without causing focal ischemia.

Built for multi-purpose intraoral manipulation, this Bishop Lip And Cheek Retractor model features an overall length of 13.5 cm (\$24) with 27 mm wide working ends. Its double-ended geometry incorporates slightly differing blade curvatures on opposing ends, allowing surgeons to adapt quickly to varying facial anatomy and surgical angles without swapping instruments. The contoured central handle provides surgical assistants with steady, non-slip manual control.

Forged from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with an anti-glare satin matte finish, Retractor Maker manufactures the Bishop Lip And Cheek Retractor to withstand rigorous hospital autoclaving cycles, delivering long-lasting clinical durability, high corrosion resistance, and ergonomic balance for dental and oral surgery suites worldwide.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractor maker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractor maker.com
retractor manufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



LinkedIn



Facebook



Instagram



X (Twitter)



YouTube



Pinterest