



RETRACTOR MAKER

PRODUCT BROCHURE OF

Castaneda Sternal Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Castaneda Sternal Retractor

SKU RM-CSR-1

\$374 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
Newborn Profile	RM-CSR-1-1	\$374
Infant Profile	RM-CSR-1-2	\$384
Adult Profile	RM-CSR-1-3	\$389

PRODUCT DESCRIPTION

The Castaneda Sternal Retractor is an internationally acclaimed cardiothoracic surgical retractor engineered specifically to provide steady, controlled sternal separation during congenital heart surgery and open cardiac reconstructions. Designed with low-profile, non-magnetic spreading arms, this instrument establishes an expansive operative window into the mediastinum and pericardial cavity while minimizing mechanical stress on fragile sternal margins.

Built for precise patient-specific sizing, the Castaneda Sternal Retractor is offered in three dedicated clinical configurations: a Newborn Profile (65.0 × 55.0 × 30.0 × 16.0 mm, \$374), an Infant Profile (72.0 × 80.0 × 45.0 × 12.0 mm, \$384), and an Adult Profile (100.0 × 125.0 × 60.0 × 15.0 mm, \$389). Each profile features anatomically contoured blade depths and curved lips that seat securely against the divided sternum. Its fine-pitch rack-and-pinion drive allows millimeter-by-millimeter expansion, preventing cartilage micro-fractures in delicate neonatal rib cages.

Manufactured from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with a reflection-free satin matte finish, Retractor Maker produces the Castaneda Sternal Retractor to withstand rigorous hospital autoclaving cycles, providing pediatric cardiac and cardiovascular surgical teams with unyielding reliability and life-saving precision.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)