



RETRACTOR MAKER

PRODUCT BROCHURE OF

Biggs Mammoplasty Retractor standard

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Biggs Mammaplasty Retractor standard

SKU RM-BMRS-1

\$69 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
Biggs Mammaplasty Retractor standard	8" (20.3 cm) length, 50.0 mm wide blade	\$69
Protective Sleeve / Box	Individual transit packaging	\$69

PRODUCT DESCRIPTION

The Biggs Mammaplasty Retractor standard is an indispensable handheld surgical retractor engineered specifically to deliver broad, atraumatic soft-tissue exposure during aesthetic and reconstructive breast surgery. Designed with an anatomically curved blade that conforms smoothly beneath parenchymal tissue flaps, this instrument enables surgical teams to maintain expansive pocket elevation while minimizing assistant fatigue.

Built for optimal tissue clearance, this Biggs Mammaplasty Retractor standard configuration features an 8-inch (20.3 cm) overall length (\$69) and a 50.0 mm widened blade profile. The wide contact face distributes mechanical counter-traction evenly across glandular and muscular planes, preventing localized tissue bruising or tearing during subpectoral pocket creation. Its ergonomic flat handle provides steady, fatigue-free manual leverage throughout prolonged dissections.

Forged from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with a reflection-free satin matte finish, Retractor Maker manufactures the Biggs Mammaplasty Retractor standard to withstand repetitive hospital sterilization cycles, providing plastic and breast reconstructive surgeons worldwide with high durability and surgical precision.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)