



RETRACTOR MAKER

PRODUCT BROCHURE OF

Tebbetts Fiber Optic Breast Augmentation Retractor,
Double Handle

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021

✓ ISO 9001:2015


✓ EN ISO 7153-1:2016

✓ cGMP Approved

✓ SCCI Registered since 2023

✓ SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Tebbetts Fiber Optic Breast Augmentation Retractor, Double Handle

SKU RM-TFOBARDH-1

\$259 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
Without Light Source — 13 selectable sizes (With/Without Teeth) — Universal fiber-optic cable adapter (Storz/Wolf/Olympus/ACMI)	RM-TFOBARDH-1-1	\$259
With Wireless LED Source — 13 selectable sizes (With/Without Teeth) — Wireless battery-operated LED light unit (AA battery, gas-sterilizable)	RM-TFOBARDH-1-2	\$335

PRODUCT DESCRIPTION

The Tebbetts Fiber Optic Breast Augmentation Retractor, Double Handle is an advanced illuminated aesthetic surgical retractor engineered specifically to provide superior manual counter-traction and shadowless optical illumination during complex breast reconstructive procedures. Designed with an ergonomic double-handle configuration, this instrument allows operating assistants to apply sustained two-handed pulling leverage across dense pectoral muscles, eliminating single-hand fatigue during lengthy pocket dissections.

Built to minimize user fatigue during long operating room sessions, this Tebbetts Fiber Optic Breast Augmentation Retractor, Double Handle model features a highly customizable sizing matrix across 13 unique blade dimensions: 80×16, 90×27, 110×25, 120×25, 130×16, 130×25, 150×25, 150×30, 150×36, 180×16, 180×25, 190×36, and 220×27 mm. Handcrafted from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic), surgeons can select between a With Teeth blade layout to prevent slippage on fibrous tissue or a Without Teeth (Smooth End) layout for atraumatic manipulation. The light source packaging offers two options: Without Light Source (\$259) for existing fiber-optic cables, or With Wireless LED Reusable Light Source (\$335), which is gas-sterilizable, runs on AA batteries, and delivers an excellent white bright light.

Forged with an anti-glare satin matte finish, Retractor Maker manufactures the Tebbetts Fiber Optic Breast Augmentation Retractor, Double Handle to withstand repetitive hospital sterilization cycles, delivering unmatched optical clarity, mechanical rigidity, and complete autoclave durability for plastic and cosmetic surgery centers worldwide.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)