



RETRACTOR MAKER

PRODUCT BROCHURE OF

Adson Self-Retaining Cerebellar Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Adson Self-Retaining Cerebellar Retractor

SKU RM-ASCR-1

\$199 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
Small — 6.5 in, 4x4 prongs	RM-ASCR-1-1	\$189
Standard — 7.5 in, 4x4 prongs	RM-ASCR-1-2	\$199
Large — 9.5 in, 4x4 prongs	RM-ASCR-1-3	\$214

PRODUCT DESCRIPTION

The Adson Self-Retaining Cerebellar Retractor is a versatile, multi-configuration neurosurgical retractor engineered to maintain stable, hands-free operative access over the posterior cranial fossa and spine. By holding paraspinal and suboccipital musculature apart under mechanical counter-traction, this instrument provides an expansive, stationary visual field while preventing assistant fatigue.

Designed to adapt precisely to varying surgical approaches, the Adson Self-Retaining Cerebellar Retractor features a uniform 19 cm body length paired with custom-selected arm and prong configurations. Surgeons can select from straight arms for shallow fields, angled arms for lateral clearance, or deep 80-degree arms, paired with either sharp prongs for firm fascial purchase or blunt prongs to protect delicate tissue margins. Its heavy-duty finger-ring ratchet mechanism maintains secure lock-down without slipping.

Forged from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with a reflection-free satin finish, Retractor Maker manufactures the Adson Self-Retaining Cerebellar Retractor to exacting hospital specifications, delivering superior tactile balance, smooth ratchet control, and complete autoclave durability.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractor maker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractor maker.com
retractor manufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)