



RETRACTOR MAKER

PRODUCT BROCHURE OF

Flexible Neck Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Flexible Neck Retractor

SKU RM-FNR-1

\$29 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
2-Prong Setup — Sharp or Blunt Ends	RM-FNR-1-1	\$29
3-Prong Setup — Sharp or Blunt Ends	RM-FNR-1-2	\$34
4-Prong Setup — Sharp or Blunt Ends	RM-FNR-1-3	\$39

PRODUCT DESCRIPTION

The Flexible Neck Retractor is an adaptable handheld surgical instrument engineered specifically to provide customized soft-tissue exposure and wound margin retraction across complex cervical and facial contours. Designed with a semi-malleable neck segment, this instrument allows surgeons and assistants to bend and contour the shaft intraoperatively to bypass visual obstructions and avoid crowding surgical instruments in tight corridors.

Built for diverse anatomical zones, the Flexible Neck Retractor features a standardized overall tool length of 16 cm that delivers balanced manual control and refined tactile feedback. The instrument is manufactured across three selectable prong configurations: a 2-prong setup (\$29), a 3-prong setup (\$34), and an expansive 4-prong setup (\$39). Each prong density is selectable with either Sharp tips for firm purchase on dense fascial layers or Blunt tips to safely protect neurovascular structures and delicate skin margins.

Manufactured from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with a reflection-free satin matte finish, Retractor Maker produces the Flexible Neck Retractor to withstand repetitive bending and standard hospital autoclave cycles, providing ENT, plastic, and general surgical teams worldwide with exceptional versatility, corrosion resistance, and complete clinical durability.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)