



RETRACTOR MAKER

PRODUCT BROCHURE OF

Leriche Retractor-Atraumatic Wide Blade

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.

Donald Cameron



Leriche Retractor-Atraumatic Wide Blade

SKU RM-LRAW-1

\$94 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
25.5 cm	RM-LRAW-1-1	\$94
and 27 cm	RM-LRAW-1-2	\$99

PRODUCT DESCRIPTION

The Leriche Retractor-Atraumatic Wide Blade is a premier deep-cavity surgical instrument engineered specifically to provide broad, non-traumatic soft-tissue and visceral displacement during major abdominal, pelvic, and thoracic procedures. Designed with an expansive, smoothly curved spatula blade, this retractor distributes manual counter-traction evenly across abdominal viscera and thick muscular walls without causing localized ischemic necrosis.

Built for dependable deep-cavity clearance, this Leriche Retractor-Atraumatic Wide Blade is handcrafted from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) and can be customized across two overall lengths: 25.5 cm (\$94) and 27 cm (\$99). Surgeons can combine these reaches with two distinct working end widths: a 45 mm wide spatula blade for standard laparotomies or an expansive 60 mm wide blade for broad visceral holding. The ergonomic flat handle provides surgical assistants with steady manual leverage throughout prolonged open surgical cases.

Forged with an anti-glare satin matte finish, Retractor Maker produces the Leriche Retractor-Atraumatic Wide Blade to withstand heavy mechanical loads and repeated hospital autoclave cycles, delivering exceptional durability, corrosion resistance, and complete autoclave durability for general and gynecological surgery centers worldwide.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)