



RETRACTOR MAKER

PRODUCT BROCHURE OF

Caspar Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Caspar Retractor

SKU RM-CR-21-1

\$64 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
Caspar Retractor	21 cm length, 17.0 × 67.0 mm blade, finger ring handle	\$64
Protective Sleeve / Box	Individual transit packaging	\$64

PRODUCT DESCRIPTION

The Caspar Retractor is a specialized handheld surgical instrument engineered specifically to provide steady, deep soft-tissue retraction and wound margin exposure during orthopedic, spinal, and reconstructive operations. Designed with an elongated right-angled blade, this instrument allows surgical assistants to maintain continuous counter-traction on dense muscle flaps with minimal manual exertion.

Built for precision deep-cavity clearance, this standard Caspar Retractor model features an overall length of 21 cm (\$64) and a 17.0 mm wide by 67.0 mm long blade. Its blade incorporates a gentle distal lip that holds fascia and muscle layers securely, preventing accidental slippage during active bone resection. The built-in ergonomic finger ring handle maximizes grip security, enhances manual control in tight clinical pathways, and significantly reduces user hand fatigue during long operative exposures.

Manufactured from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with a reflection-free satin matte finish, Retractor Maker produces the Caspar Retractor to strict clinical standards, delivering high tensile strength, exceptional corrosion resistance, and complete autoclave durability for hospital operating theaters worldwide.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)