



RETRACTOR MAKER

PRODUCT BROCHURE OF

Neonatal T Ribbon Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Neonatal T Ribbon Retractor

SKU RM-NTRR-265-1

\$49 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
2.5 cm × 1.25 cm	RM-NTRR-265-1-1	\$49
3.8 cm × 1.40 cm	RM-NTRR-265-1-2	\$54
4.5 cm × 1.60 cm	RM-NTRR-265-1-3	\$59

PRODUCT DESCRIPTION

The Neonatal T Ribbon Retractor is a specialized pediatric and neonatal surgical instrument engineered specifically to provide gentle, atraumatic visceral displacement and wound margin exposure in fragile infant operating fields. Designed with a flexible malleable shaft, this retractor allows pediatric surgeons to contour the blade intraoperatively to match delicate neonatal abdominal and thoracic curves without exerting localized pressure points.

Built with a uniform total reach framework of 26.5 cm, this Neonatal T Ribbon Retractor accommodates varying surgical fields through three selectable blade footprint dimensions: 2.5 cm × 1.25 cm (\$49), 3.8 cm × 1.40 cm (\$54), and 4.5 cm × 1.60 cm (\$59). Handcrafted from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic), its working blade geometry features a malleable, T-shaped end configuration profile with rounded, atraumatic margins that protect infant viscera from contusions. The slender shank provides surgical assistants with fatigue-free manual leverage.

Forged with an anti-glare satin matte finish, Retractor Maker manufactures the Neonatal T Ribbon Retractor to withstand repeated intraoperative re-contouring and hospital sterilization cycles, providing pediatric surgery centers and congenital heart institutes worldwide with exceptional ductility, corrosion resistance, and complete autoclave durability.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)