



RETRACTOR MAKER

PRODUCT BROCHURE OF

Ollier Retractors

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Ollier Retractors

SKU RM-OR-SERIES-1

\$74 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
8 ½" (21.59 cm)	RM-OR-SERIES-1-1	\$74
9.0" (22.86 cm)	RM-OR-SERIES-1-2	\$79

PRODUCT DESCRIPTION

Ollier Retractors are foundational heavy-duty surgical instruments engineered specifically to provide broad, atraumatic soft-tissue displacement and muscular retraction during major open abdominal, pelvic, and orthopedic procedures. Designed with an expansive 4-prong blunt rake head, these retractors allow surgical assistants to apply heavy manual counter-traction across thick muscular walls and subcutaneous fat without causing focal tissue contusions or lacerations.

Built for lasting leverage in deep operative spaces, the Ollier Retractors line keeps a uniform multi-pronged rake profile while offering two selectable lengths from the purchasing menu parameters: 8 ½ inches (21.59 cm, \$74) and 9 inches (22.86 cm, \$79). Handcrafted from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic), the structural geometry features a 4 blunt pronged rake configuration coupled with a hollow comfort handle that significantly reduces tool weight while maximizing grip stability, preventing assistant hand and wrist fatigue during multi-hour procedures.

Forged with an anti-glare satin matte finish, Retractor Maker manufactures Ollier Retractors to withstand heavy mechanical loads and repeated hospital autoclave cycles, providing general, orthopedic, and trauma surgical teams worldwide with exceptional durability, corrosion resistance, and complete clinical dependability.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)