



RETRACTOR MAKER

PRODUCT BROCHURE OF

Biggs Mammoplasty Retractor narrow

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Biggs Mammaplasty Retractor narrow

SKU RM-BMRN-1

\$89 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
Biggs Mammaplasty Retractor narrow	6" (15.2 cm) length, 3 cm narrow curved blade	\$89
Protective Sleeve / Box	Individual transit packaging	\$89

PRODUCT DESCRIPTION

The Biggs Mammaplasty Retractor narrow is an essential handheld plastic surgery retractor engineered specifically to provide clear, atraumatic soft-tissue exposure during aesthetic and reconstructive breast surgery. Designed with an anatomical profile, this instrument allows surgeons to elevate parenchymal and muscular flaps with optimal leverage while minimizing counter-traction strain.

Built for precision breast pocket dissection, this Biggs Mammaplasty Retractor narrow model features an overall length of 6 inches (15.2 cm) and a compact 3 cm wide blade profile. The strongly curved shaft is engineered to adapt smoothly around anatomical breast curves, giving surgeons a clear, unobstructed view of deep subpectoral and subglandular pockets without causing localized soft-tissue bruising or tearing.

Forged from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with a reflection-free satin matte finish, Retractor Maker manufactures the Biggs Mammaplasty Retractor narrow to withstand repetitive hospital sterilization cycles, providing breast and plastic surgeons worldwide with superior ergonomic balance, corrosion resistance, and complete autoclave durability.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)