



RETRACTOR MAKER

PRODUCT BROCHURE OF

Johnson Neuroma Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Johnson Neuroma Retractor

SKU RM-JNR-1

\$74 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
8.0 cm	RM-JNR-1-1	\$74
10.0 cm	RM-JNR-1-2	\$79
11.0 cm	RM-JNR-1-3	\$84
14.0 cm	RM-JNR-1-4	\$89

PRODUCT DESCRIPTION

The Johnson Neuroma Retractor is a specialized handheld surgical instrument engineered specifically to provide gentle, non-traumatic soft-tissue and metatarsal head retraction during Morton's neuroma excision and interdigital nerve procedures. Designed to work in confined web-space corridors, this instrument allows podiatric and orthopedic surgeons to mobilize intermetatarsal ligaments and digital vessels without crushing fragile nerve trunks.

Built with highly adaptive size ranges and profile varieties, the Johnson Neuroma Retractor is available across four overall lengths: 8.0 cm (\$74), 10.0 cm (\$79), 11.0 cm (\$84), and 14.0 cm (\$89). Surgeons can choose between three design typologies: a Narrow profile layout for tight intermetatarsal spaces, a 90° Angled profile with a 7.5 mm wide working blade, or a 90° Angled profile with an expanded 9.5 mm wide working blade. The ergonomic flat handle provides tactile sensitivity and non-slip control during delicate digital neurolysis.

Handcrafted from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with a reflection-free satin matte finish, Retractor Maker produces the Johnson Neuroma Retractor to strict clinical tolerances, delivering absolute corrosion resistance, unyielding structural rigidity, and complete autoclave durability for foot and ankle surgery suites worldwide.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)