



RETRACTOR MAKER

PRODUCT BROCHURE OF

Salzburg Mandibular Rim Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Salzburg Mandibular Rim Retractor

SKU RM-SMRR-23-1

\$49 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
Salzburg Mandibular Rim Retractor	23 cm overall length, curved mandibular rim blade	\$49
Protective Sleeve / Box	Individual transit packaging	\$49

PRODUCT DESCRIPTION

The Salzburg Mandibular Rim Retractor is an indispensable maxillofacial surgical retractor engineered specifically to provide safe, non-traumatic soft-tissue displacement and osseous border exposure along the lower border of the mandible during orthognathic interventions. Designed with an anatomical rim curve that hooks smoothly around the cortical bone of the mandibular corpus and angle, this instrument allows oral surgeons to perform osteotomy bone cuts cleanly while shielding adjacent neurovascular structures from reciprocating saws.

Built for optimal intraoral reaching clearance, this Salzburg Mandibular Rim Retractor configuration features an overall tool length of 23 cm (\$49). Handcrafted from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) to ensure lasting strength and corrosion resistance, its smoothly polished blade margins cradle the lower mandibular border without causing localized periosteal crush trauma or lacerating deep cervical tissues. The ergonomic flat handle provides operating assistants with steady, fatigue-free manual counter-traction throughout extended orthognathic procedures.

Forged with an anti-glare satin matte finish, Retractor Maker manufactures the Salzburg Mandibular Rim Retractor to withstand repetitive hospital autoclave cycles, delivering exceptional mechanical stability, corrosion resistance, and complete autoclave durability for oral, maxillofacial, and plastic surgery suites worldwide.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)