



RETRACTOR MAKER

PRODUCT BROCHURE OF

Lobster Anterior Lumbar Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.

Donald Cameron



Lobster Anterior Lumbar Retractor

SKU RM-LALR-1

\$3,999 USD

COMPONENT LIST

QUANTITY	COMPONENT	DIMENSIONS / DESCRIPTION
6	Lobster Retractor Clamp	—
1	Lobster Table Clamp with Vertical Post & Frame Coupler	—
1	Lobster Frame	—
2	Lobster Mayo Swivel Retractor, Stainless Steel	2 ³ / ₄ " × 2"
1	Lobster Mayo Swivel Retractor, Stainless Steel	3 ¹ / ₂ " × 2"
1	Lobster Deaver Swivel Retractor, Stainless Steel	2 ³ / ₄ " × 4 ¹ / ₂ "
1	Lobster Harrington Swivel Retractor, Stainless Steel	2" × 5"
1	Lobster Malleable Swivel Retractor, Stainless Steel	2" × 6 ¹ / ₂ "
1	Lobster Malleable Swivel Retractor, Stainless Steel	3" × 6 ¹ / ₂ "
1	Lobster Fence Swivel Retractor, Stainless Steel	4" × 6"

QUANTITY	COMPONENT	DIMENSIONS / DESCRIPTION
2	Lobster Renal Swivel Retractor, Stainless Steel	1" × 4"
2	Lobster Renal Swivel Retractor, Stainless Steel	1" × 6"
2	Lobster Renal Swivel Retractor, Stainless Steel	1" × 8"
1	Lobster Splanchnic Swivel Retractor, Stainless Steel	1½" × 6"
1	Lobster Splanchnic Swivel Retractor, Stainless Steel	2" × 4½"
1	Lobster Splanchnic Swivel Retractor, Stainless Steel	2" × 6"
1	Lobster Splanchnic Swivel Retractor, Stainless Steel	2" × 8"
1	Sterilization Box	—

PRODUCT DESCRIPTION

The Lobster Anterior Lumbar Retractor is a heavy-duty table-mounted retractor system designed specifically to address the demanding exposure challenges of anterior approach spine surgery. Mounting directly to the operating table accessory rail via a solid post and frame coupler, this rigid platform eliminates assistant fatigue while providing stationary, reliable multi-planar retraction across the retroperitoneal space.

Engineered to protect critical vascular structures during disc space preparation, the Lobster Anterior Lumbar Retractor system incorporates a rigid Lobster structural frame and six quick-locking retractor clamps. The set comes furnished with an extensive 16-piece swivel retractor blade complement—including narrow Renal blades (up to 8" deep), broad Splanchnic visceral retractors, protective Fence blades, Deaver, Harrington, Mayo, and bendable Malleable blades that swivel smoothly to adapt to patient-specific lordosis and anatomy.

Manufactured from biocompatible, Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with an anti-reflective satin finish, Retractor Maker supplies the Lobster Anterior Lumbar Retractor complete in a heavy-duty sterilization box, providing hospitals and spinal centers with a fully autoclavable, high-precision solution for complex lumbar reconstructions.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)