



RETRACTOR MAKER

PRODUCT BROCHURE OF

Aufricht Daniel Nasal Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Aufricht Daniel Nasal Retractor

SKU RM-ADNR-1

\$29 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
Aufricht Daniel Nasal Retractor	11 cm overall length, 50 × 10 mm blade	\$29
Protective Sleeve / Box	Individual transit packaging	\$29

PRODUCT DESCRIPTION

The Aufricht Daniel Nasal Retractor is a specialized, compact rhinoplasty retractor engineered to provide clear, atraumatic visualization of the nasal dorsum, cartilaginous vault, and septal framework. Designed with a shortened 11 cm profile, this instrument allows operating surgeons and assistants to achieve high maneuverability and steady tissue elevation in confined operative spaces.

Optimized to reduce intraoperative tissue strain, this Aufricht Daniel Nasal Retractor configuration features a 50 mm long by 10 mm wide smooth blade that slips easily under the dorsal skin envelope. The blade assembly incorporates a gentle downward angle that conforms to the nasal bridge, providing broad-surface contact that eliminates localized pressure points. Its proportionate, lightweight handle reduces hand fatigue during lengthy aesthetic nasal remodeling.

Forged from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with a reflection-free satin matte finish, Retractor Maker manufactures the Aufricht Daniel Nasal Retractor to strict clinical standards, ensuring superior corrosion resistance, easy reprocessing, and complete autoclave compatibility for plastic and ENT surgical suites.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)