



RETRACTOR MAKER

PRODUCT BROCHURE OF

Cerkes Double Ended Nasal Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Cerkes Double Ended Nasal Retractor

SKU RM-CDENR-1

\$69 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
Small (10×7 mm & 15×7 mm)	RM-CDENR-1-1	\$69
Large (30×8 mm & 35×10 mm)	RM-CDENR-1-2	\$74

PRODUCT DESCRIPTION

The Cerkes Double Ended Nasal Retractor is a specialized rhinoplasty surgical retractor engineered specifically to provide multi-depth tissue elevation and precise anatomical exposure during cosmetic and reconstructive nasal procedures. Designed in collaboration with advanced rhinoplasty techniques, this dual-ended instrument enables surgeons to mobilize the nasal skin envelope and lower lateral cartilages without switching instruments mid-dissection.

Built for refined procedural flexibility, this Cerkes Double Ended Nasal Retractor model features an overall length of 15 cm and is manufactured in two selectable dual-ended blade configurations. The Small Configuration (\$69) integrates a 10 mm × 7 mm blade on one side and a 15 mm × 7 mm blade on the opposite side, while the Large Configuration (\$74) features an expansive 30 mm × 8 mm blade paired with a 35 mm × 10 mm blade. The smooth, contoured blade edges prevent mucosal tears while maintaining rigid, stable retraction.

Forged from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with an anti-glare satin matte finish, Retractor Maker produces the Cerkes Double Ended Nasal Retractor to withstand repetitive hospital autoclave cycles, providing plastic and ENT surgeons worldwide with exceptional tactile balance, corrosion resistance, and surgical longevity.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)