



RETRACTOR MAKER

PRODUCT BROCHURE OF

Femoral/Gluteus Medius Minimus Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Femoral/Gluteus Medius Minimus Retractor

SKU RM-FGR-1175-1

\$44 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
Femoral/Gluteus Medius Minimus Retractor	11.75" (29.85 cm) length, 20 mm wide blade	\$44
Protective Sleeve / Box	Individual transit packaging	\$44

PRODUCT DESCRIPTION

The Femoral/Gluteus Medius Minimus Retractor is a heavy-duty orthopedic surgical instrument engineered specifically to handle dense, bulky muscular structures during total hip arthroplasty and complex pelvic reconstruction. Designed with an extended reach and an anatomically contoured blade, this retractor allows surgical assistants to leverage the gluteus medius, gluteus minimus, and anterior capsular tissues away from the joint corridor without causing muscle tears.

Built to handle heavy muscle structures with ease, this Femoral/Gluteus Medius Minimus Retractor configuration features an overall tool length of 11.75 inches (29.85 cm, \$44) paired with a 20 mm wide blade profile. Its gently curved blade geometry wraps smoothly around the superior and posterior aspects of the greater trochanter and femoral neck, preventing accidental slippage off bone cortices during deep acetabular reaming. The ergonomic handle provides superior pulling leverage, reducing assistant hand strain.

Forged from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with a reflection-free satin matte finish, Retractor Maker produces the Femoral/Gluteus Medius Minimus Retractor to withstand intense mechanical torque and repeated hospital autoclave cycles, providing orthopedic surgery teams worldwide with exceptional durability and clinical performance.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)