



RETRACTOR MAKER

PRODUCT BROCHURE OF

Austin Tissue and Cheek Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Austin Tissue and Cheek Retractor

SKU RM-ATCR-1

\$45 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
Austin Tissue and Cheek Retractor	Compact right-angled intraoral retractor	\$45
Protective Sleeve / Box	Individual transit packaging	\$45

PRODUCT DESCRIPTION

The Austin Tissue and Cheek Retractor is a specialized intraoral surgical instrument engineered to provide steady, non-slip retraction of the buccal mucosa, lips, and mucoperiosteal flaps during dental and oral surgical procedures. Designed with an ergonomic right-angled profile, this retractor securely holds soft tissues away from surgical drills and rotary burs without causing tissue ischemia.

Built for optimal intraoral clearance, this Austin Tissue and Cheek Retractor model features a flat, smooth right-angled blade tailored to seat comfortably within the oral vestibule. Its low-profile architecture ensures that the retractor blade does not obstruct the surgeon's operative view or impinge on dental arches. The balanced handle provides surgical assistants with fatigue-free manual control during prolonged surgical extractions and implant placements.

Forged from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with a reflection-free satin matte finish, Retractor Maker manufactures the Austin Tissue and Cheek Retractor to strict clinical standards, ensuring superior corrosion resistance, easy cleaning, and complete autoclave compatibility for oral surgery and dental clinics worldwide.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)